

TEV protease - (His-Strep)N-tag, endotoxin-free

1.0 mg, 1.01 mg/ml

Catalog UCT004

Lot 2708

product specification sheet: T004



products@immunoprecise.com

www.immunoprecise.com

Product description

The Tobacco etch virus (TEV) protease is one of the most popular proteases to remove affinity tags from a recombinant protein containing the TEV-protease recognition site Glu-Asn-Leu-Tyr-Phe-Gln-Gly. Cleavage occurs between the Gln and Gly residues, resulting in a Gly at the amino terminus of the protein of interest.

The TEV protease is active in 2.5% sucrose, 0.01% TRITON X-100, PMSF, aprotinin, and leu-peptin. However, 0.01 % SDS completely inhibits TEV.

Reference sequences: NC001555
NP_062908
NP_734212
Swiss-prot P04517

Alternative names: Peptidase_C4
Nla-Pro protein
Tobacco etch virus protease

TEV protease is produced with an N-terminal (His-Strep)-tag and purified to homogeneity. It is produced without any animal components. The N-terminal (His-Strep)-tag facilitates efficient removal of the protease by affinity chromatography after digestion is completed.

Protocol

For each application most optimal digestion conditions should be determined experimentally. For large scale digestions a physiological buffer should be used, for example PBS.

To start with use the following digestion conditions:

- * TEV-(Strep-His)N protease: 0.02 mg/ml working concentration
- * Substrate (10-100 kDa): 1 - 5 mg/ml
- * Incubation temperature: room temperature
- * Incubation time : 1 hour - over night
- * Digestion buffer: PBS

Digestion should be performed with gentle mixing.

Storage, stability and characteristics

Recombinant TEV-(His-Strep) protease is sterile and should be stored at -20 °C (stable for at least 6 months). The buffer contains PBS, 1 mM DTT, 50% glycerol without a preservative.

Endotoxin: < 0.4 EU/mg as determined by a LAL assay

Purity: > 98 %

PRODUCT USE LIMITATIONS, WARRANTY, DISCLAIMER

ImmunoPrecise Antibodies Europe BV products contain chemicals which may be harmful if misused. Due care should be exercised with all ImmunoPrecise Antibodies Europe BV products to prevent direct human contact. All products are intended For Research Use Only and ARE NOT ALLOWED FOR USE IN HUMANS. Each ImmunoPrecise Antibodies Europe BV product is shipped with documentation stating specifications and other technical information. ImmunoPrecise Antibodies Europe BV products are warranted to meet or exceed the stated specifications. ImmunoPrecise Antibodies Europe BV's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. ImmunoPrecise Antibodies Europe BV makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO ImmunoPrecise Antibodies Europe BV PRODUCTS. In no event shall ImmunoPrecise Antibodies Europe BV be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of ImmunoPrecise Antibodies Europe BV products to perform in accordance with the stated specifications.