

Recombinant FSH (Danio rerio, single chain)

Catalog F001
100 µg
Lot A00017



H.R. Kruytgebouw
Padualaan 8
3584 CH Utrecht
The Netherlands

Product description

The pituitary glycoprotein hormone family includes follicle-stimulating hormone, luteinizing hormone, chorionic gonadotropin, and thyroid-stimulating hormone. All of these glycoproteins consist of an identical alpha subunit and a hormone-specific beta subunit that confers biological specificity.

Reference sequences: NM_205624.1 and NM_205687
NP_991187 and NP_991250
Swiss-prot Q6TCF7 and Q6QX43

Alternative names: follitropin subunit beta precursor
follitropin subunit alpha precursor

Gene ID: 402919 and 402987

This protein is encoded by a fusion of the mature follitropin subunit beta (residues 17 - 130) and glycoprotein hormones alpha chain (residues 24 - 117) linked by a (gggs)₄ peptide (Nature, 2005, vol 433, 269-277). It contains a N-terminal Strep-tag. The protein was over-expressed in HEK293EBNA1 cells and purified to homogeneity (figure 1). The calculated molecular weight of recombinant Danio rerio FSH is 27.9 kDa. Each vial contains 100 µg FSH at a concentration of 0.98 mg/ml.

Protein sequence

gswshpqfekgswshpqfekgswshpqfekgsaesecrcsrltnisitveseecgscvtidttacaglcwtmdrvypssmaqhtqkvcnfknImyksyefkgcpagvds
vfvyvpalscecnqvnsttdtdwgaispqttsctsihgsgsgsgsgsgsgsgsrndvsnygceeklkmmnerfskpgapvyqcvgccfsrayptplrskktmlvpknitse
atccvakeskmvatniplynhtdchcstcyhks

The Strep-tag is shown in red.

Storage and stability

FSH should be stored at - 80 °C (stable for at least 1 year). The buffer contains PBS without preservative. After thawing it should be stored in appropriate small aliquots at - 20 °C or - 80 °C (stable for at least 2 months).




Figure 1.
NuPage analysis of purified FSH-StrepN.
lane 1: marker
lane 2: FSH-StrepN

PRODUCT USE LIMITATIONS, WARRANTY, DISCLAIMER

U-Protein Express BV products contain chemicals which may be harmful if misused. Due care should be exercised with all U-Protein Express BV products to prevent direct human contact. All products are intended For Research Use Only and ARE NOT ALLOWED FOR USE IN HUMANS. Each U-Protein Express BV product is shipped with documentation stating specifications and other technical information. U-Protein Express BV products are warranted to meet or exceed the stated specifications. U-Protein Express BV's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. U-Protein Express BV makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO U-Protein Express BV PRODUCTS. In no event shall U-Protein Express BV be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of U-Protein Express BV products to perform in accordance with the stated specifications.